

PEDUS

Fore- and Midfoot Plating System

PEDUS MTP 1 Fusion Plate

Indicated for fixation of fractures, osteotomies and joint fusions at the fore- and midfoot, especially at the first metatarsophalangeal joint with Hallux Rigidus.

PEDUS-L Locking Lapidus Plate

Indicated for the fixation of fractures, osteotomies and joint fusions at the midfoot, especially at the first tarsometatarsal joint. The plates may be used for correction of deformities, especially hallux valgus.

PEDUS-Lx Compression Plate

Indicated for the fixation of fractures, osteotomies and joint fusions at the midfoot, especially at the first tarsometatarsal joint. The plates may be used for correction of deformities, especially hallux valgus. The lag screw enables secure screw fixation into the Os cuneiforme 1 or 2.

PEDUS-MTP1

Fusion Plate

Laser marking

For alignment over the metatarsophalangeal joint to be fused

Screw holes suitable for:


- Cortical Screws
- Locking Cortical Screws
- Multiaxial Locking Cortical Screws

Anatomically shaped
With 5° valgus angle

Dorsal extension 0° or 8

Low Profile

Protection of the soft tissue environment



PEDUS-Lx

Compression Plate

Laser marking

For alignment over the TMT1 joint to be fused

Screw holes suitable for:

- Cortical Screws
- Locking Cortical Screws
- Multiaxial Locking Cortical Screws

5° and 9° valgus angles

Lag screw

Deflection in Os cuneiforme 1 or 2 possible



PEDUS-L

Plantar Lapidus Plate

Two different step heights
3 und 5 mm

Screw holes suitable for:

- Cortical Screws
- Locking Cortical Screws
- Multiaxial Locking Cortical Screws





PEDUS-L Locking Lapidus Plate

12.11123.XXX(S)

Holes 4

Steps 0, 2, 3, 4, 5 mm

Length 31 mm



PEDUS-L Locking Plantar Lapidus Plate

12.11124.XXX(S)

Holes 4

Steps 3, 5 mm

Length 39, 45 mm



PEDUS Locking T-Plate

12.11210.XXX(S)

Holes Plate Head 2

Holes Plate Shaft 3, 4, 6

Length 32, 39, 54 mm



PEDUS Locking Plate, straight

12.11211.XXX(S)

Holes 4, 6, 8

Length 30, 43, 56 mm



PEDUS Locking L-Plate, 90°

12.11212.1XX(S) left

12.11212.0XX(S) right

Length Plate Head 2

Length Plate Shaft 3, 4, 6

Length 32, 40, 55 mm



PEDUS-Lx Compression Plate

12.11230.1XXS left

12.11230.0XXS right

Extension 5°, 9°

Length 37.50 mm



PEDUS MTP 1 Fusion Plate, short

12.1125X.303S left

12.1125X.203S right

Holes Plate Head 2

Holes Plate Shaft 3

Extension 0°, 8°

Length 37 mm



PEDUS MTP 1 Fusion Plate

12.1125X.10X(S) left

12.1125X.00X(S) right

Holes Plate Head 3

Holes Plate Shaft 3, 5

Extension 0°, 8°

Length 53, 69 mm



PEDUS Locking Multifix Plate

12.11018.007(S)

Holes 7

Diameter 18 mm



PEDUS Locking X-Plate


12.11220.XXX(S)

Holes 4


Length 27, 33 mm

PEDUS

Fore- and Midfoot Plate System



PEDUS Instruments



Cortical Screw Ø 2.7 mm, self-tapping

03.03527.xxx(S)

Length 08-40 mm



Locking Cortical Screw Ø 3.0 mm, self-tapping

10.03530.xxx(S)

Length 08-40 mm



Washer Ø 7.0 mm, for Screws Ø 2.7 - 4.0 mm

03.91000.070

Outer diameter 7.0 mm



PEDUS Templates

12.11xxx.xxx

1 = Implant

12.21xxx.xxx

2 = Template

Templates are available for all plate types



Multiaxial Locking Cortical Screw Ø 2.7 mm, self-tapping

03.03540.xxx(S)

Length 08-36 mm



Cannulated Screw Ø 4.0 mm, short thread, self-drilling

08.03644.xxx(S)

Length 20-50 mm

You can also find us on social networks:



Note:
This description is not sufficient for immediate application of the instruments. Further information is available in the detailed surgical technique.
Maintenance, care and preparation of Marquardt instruments can be found in the corresponding instructions.